



QuickGene Series Application Guide

Total RNA Extraction from Cultured Cells Total RNA Extraction from Cells Cultured in 6cm, 10cm Dish

Kit : QuickGene SP kit RNA cultured cell HC (Spin method)

Protocol


< Number of cells and corresponding protocols >

Before performing an extraction, please select an extraction protocol by reference to following table.

Dish size	Cell species	Number of cells (x 10 ⁶)	Protocol
6cm	HeLa	1.5~3.5	A
	HEK293	3.0~5.0	
	COS-7	0.5~1.5	
	NIH/3T3	1.0~2.5	
	HL60	3.0~5.0	
10cm	HeLa	3.5~5.5	B
	HEK293	5.0~8.0*	
	COS-7	2.0~3.0	
	NIH/3T3	3.0~5.0	
	HL60	5.0~15	
	HEK293	8.0~15*	B'


* If lysing a number of cells exceeding 8.0 x 10⁶ directly on a dish, use Protocol B'.
Pelleted cells can be processed up to 15 x 10⁶ with Protocol B.

● Protocol A




*1 : 2-Mercaptoethanol (2-ME) must be added to LRP before each use. Add 10 µl 2-ME per 1 ml of LRP.

● Protocol B




*1 : 2-Mercaptoethanol (2-ME) must be added to LRP before each use. Add 10 µl 2-ME per 1 ml of LRP.

● Protocol B'



*1 : 2-Mercaptoethanol (2-ME) must be added to LRP before each use. Add 10 µl 2-ME per 1 ml of LRP.



*2 : Add 25 ml of >99% ethanol into the bottle and mix by gently inverting the bottle before use.

*3 : DNase is not included in the kit. Please prepare following recommended product.

Recommended DNase

- | | |
|---------------------------------|------------------------------------|
| a) RQ1 RNase-Free DNase | (Promega : Cat. No. M6101) |
| b) DNase I, Amplification Grade | (Invitrogen : Cat. No. 18068-015) |
| c) DNase I, Amplification Grade | (SIGMA : Cat. No. AMP-D1) |
| d) Deoxyribonuclease (RT Grade) | (Nippon Gene : Cat. No. 313-03161) |
| e) DNase I, RNase-Free | (Ambion : Cat. No. 2222) |
| f) RNase-Free DNase Set | (QIAGEN : Cat. No. 79254) |

In the case of DNase a)~d)

1 U / µl DNase I	: 20 µl
10 x Reaction Buffer	: 4 µl
Nuclease-free water	: 16 µl

In the case of DNase e)

2 U / µl DNase I	: 20 µl
10 x Reaction Buffer	: 4 µl
Nuclease-free water	: 16 µl

In the case of DNase f)

2.7 Kunitz unit / µl DNase I	: 1.25 µl
Buffer RDD	: 35 µl
Nuclease-free water	: 3.75 µl

Results : Total RNA extraction from various model cells


Total RNA was extracted from various model cells with QuickGene SP kit RNA cultured cell HC.

● The yield and purity of total RNA

Cells species		HL60		HEK293		HeLa		COS-7		NIH / 3T3	
Protocol		A	B	A	B'	A	B	A	B	A	B
Cell form		Pellet		6 cm dish	10 cm dish	6 cm dish	10 cm dish	6 cm dish	10 cm dish	6 cm dish	10 cm dish
Number of cells (x 10 ⁶ cells)		5.0	15	4.5	9.9	1.7	4.9	1.0	2.7	1.7	4.5
DNase(+)	Yield(μg)	44.0	167.4	89.1	213.0	52.0	126.5	40.4	110.0	32.5	97.1
	A _{260/280}	2.17	2.15	2.18	2.16	2.17	2.19	2.21	2.21	2.23	2.20
	A _{260/230}	2.23	2.18	1.98	2.21	2.08	2.23	2.06	2.22	2.15	2.21
DNase(-)	Yield(μg)	47.7	165.7	84.0	224.2	-	127.7	-	-	-	-
	A _{260/280}	2.17	2.15	2.19	2.17	-	2.18	-	-	-	-
	A _{260/230}	2.22	2.21	2.01	2.21	-	2.13	-	-	-	-


● Electrophoresis of total RNA

Electrophoresis conditions : Non denaturing gel electrophoresis (1% Agarose gel / 1 x TAE Buffer)



● RT-PCR (with DNase treatment)

RT-PCR was performed with total RNA extracted from various sample cells using QuickGene SP kit RNA cultured cell HC.



* Trademark and exclusion item

Right to registered names etc. used in this Application Guide is protected by law especially even in the case of no denotation.

FUJIFILM