

# FUJINON Video Choledocho Endoscopes



**FUJINON**  
 FUJIFILM

## FUJINON VIDEO CHOLEDOCHO ENDOSCOPES

# EO-270F

Fujinon is renowned for its versatility and advanced technology in all areas of endoscopy. In order to meet the diverse needs of today's healthcare providers, Fujinon has developed a new specialty scope in the field of General Surgery, the EO-270F Video Choledochoscope.



- The EO-270F is a truly unique endoscope that incorporates a large 2.2mm working channel within a slim 5.1mm insertion section. The slim size increases maneuverability and allows for easy introduction even through areas with severe strictures.
- Thanks to an ultra-small, high density CCD chip, the EO-270F enables you to see the highest resolution images. This advancement results in outstandingly crisp images. In addition, the level of observation and overall diagnostic potential has been improved on the video monitor.
- In order for the endoscope control section to be held more comfortably, a newly designed lightweight control section has been designed. This innovation also lessens hand fatigue during the operation.

	EO-270F		
Viewing direction	Forward		
Observation range	3 ~ 100mm		
Field of view	120°		
Distal end diameter	5.1mm		
Flexible portion diameter	5.3mm		
Bending capability : UP	180°		
	DOWN		130°
Forceps channel diameter	2.2mm		
Working length	450mm		
Total length	720mm		

FUJIFILM GROUP

### FUJINON CORPORATION

1-324 UETAKE, KITA-KU, SAITAMA CITY, SAITAMA 331-9624, JAPAN  
 TEL : 81-48-668-2153 FAX : 81-48-668-1570  
 E-mail:sales@msv.fujinon.co.jp URL <http://www.fujinon.co.jp/>

### FUJINON INC.

10 HIGH POINT DRIVE, WAYNE, NJ 07470, U.S.A.  
 TEL : 1-800-490-0661 FAX : 1-973-633-8818  
 E-mail : med.info@fujinon.com URL <http://www.fujinonendoscopy.com/>

### FUJINON (EUROPE) GmbH

HALSKESTRASSE 4, 47877 WILLICH, GERMANY  
 TEL : 49-2154-924-0 FAX : 49-2154-924-290  
 E-mail : Fujinon@fujinon.de URL <http://www.fujinon.de/>

### FUJINON SINGAPORE PTE. LTD.

BLK 211 HENDERSON ROAD, #10-04 HENDERSON INDUSTRIAL PARK,  
 SINGAPORE, 159552 TEL : 65-6276-4988 FAX : 65-6276-6911  
 E-mail : fujinon@fujinon.com.sg URL <http://www.fujinon.com.sg/>

### FUJINON AUSTRALIA PTY. LTD.

UNIT-18, 52 HOLKER STREET, SILVERWATER N.S.W. 2128, AUSTRALIA  
 TEL : 61-2-9748-2744 FAX : 61-2-9748-2428 E-mail : sales@fujinon.com.au

### SHANGHAI FUJINON MEDICAL EQUIPMENT CO.,LTD.

FLAT B, 6/FL, SHI MEI BLD,79 RI JING ROAD,  
 WAI GAO QIAO FREE TRADE ZONE, SHANGHAI, ZIP CODE: 200131 CHINA  
 TEL : 86-21-5424-9262(Office) TEL:86-21-5866-2151(Service)  
 FAX : 86-21-5866-2150 E-mail : fjng2001@hotmail.com