


EL-450FP

Flexible Video Laparoscope

EL-450FP



Outline



Specifications

Distal end diameter	5.4mm	
Rigid portion diameter	5.4mm	
Working length	350mm	
Total length	570mm	
Bending capability	UP	100°
	DOWN	100°
	RIGHT	100°
	LEFT	100°
Viewing direction	Forward	
Field of view	80°	
Observation range	18~120mm	

● Processor VP-4400 · Light Source XL-4400



Super slim diameter
5.4mm

The unique flexible distal end of the scope provides an optimum angle of view. The flexible distal end allows the surgeon to aim at the area of interest during complicated manipulation of the instrument inside the abdominal cavity.